



NobelProcera™ Abutment Titanium and Zirconia

INSTRUCTIONS FOR USE

22671 A West 1011 Printed in Sweden. Date of issue 2009-10-13. All rights reserved. Nobel Biocare, the Nobel Biocare logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Nobel Biocare. Product images in this folder are not necessarily to scale.

	Procera® Abutment Titanium		Procera® Abutment Zirconia	
	Abutment screw	Lab screw	Abutment screw	Lab screw
Repl NP	29474	31170	28837	-
Repl RP, WP, 6.0	29475	29293	28816	29293
Brånemark NP	29282	31168	31171	31168
Brånemark RP	29283	29290	28815	29290
Brånemark WP	29284	31169	28844	31169
NobelActive™ NP	31171	31168	31171	31168
NobelActive™ RP	28815	29290	28815	29290

NobelProcera™ Abutment Titanium and Zirconia screw chart (table 2).	
Brånemark System®, NobelReplace™, NobelActive™	35 Ncm
Astra Tech® Yellow	15 Ncm
Astra Tech® Aqua	20 Ncm
Astra Tech® Lilac	25 Ncm
Straumann®	35 Ncm

NobelProcera™ Abutment Titanium and Zirconia screw chart	
Procera® Abutment Titanium	Procera® Abutment Zirconia
Abutment screw	Abutment screw
Lab screw	Lab screw

Argentina Nobel Biocare Argentina Phone: +54 11 4625 9696 Cust. support: toll free 0800 800 66235	Belgium Nobel Biocare Belgium Phone: +32 2 467 41 70 Cust. support: toll free 0800 800 66235	Brazil Nobel Biocare Brazil Phone: +55 11 5102 7000 Cust. support: 0800 169 996	Canada Nobel Biocare Canada Phone: +1 905 762 3500 Cust. support: toll free 1800 804 597	Australia Nobel Biocare Australia Phone: +61 2 8064 5100 Cust. support: toll free 1800 804 597	Austria Nobel Biocare Austria Phone: +43 1 892 89 90
--	---	---	---	---	---

China Nobel Biocare China Phone: +86 21 5206 6655 Cust. support: +86 21 5206 0974	Denmark Nobel Biocare Denmark Phone: +45 39 40 48 46	Finland Nobel Biocare Finland Phone: +358 9 343 69 70	France Nobel Biocare France Phone: +33 1 49 20 00 30 Cust. support: +86 21 5206 0974	Germany Nobel Biocare Germany Phone: +49 221 500 85 980	Hong Kong Nobel Biocare Hong Kong Phone: +852 2845 1266 Cust. support: +852 2823 8926
--	---	--	---	--	--

Hungary Nobel Biocare Hungary Phone: +36 1 279 33 79	India Nobel Biocare India Phone: +91 22 6751 9999 Cust. support: toll free 1 800 22 9998	Israel Nobel Biocare Israel Phone: +972 874 5951	Italy Nobel Biocare Italy Phone: +39 039 683 61 Cust. support: toll free 800 53 93 28	Japan Nobel Biocare Japan Phone: +81 3 6717 6191	Lithuania Nobel Biocare Lithuania Phone: +370 5 269 3448 Cust. support: toll free 800 01 23 24
---	---	---	--	---	---

Middle East Nobel Biocare Middle East Phone: +48 22 874 5951	Netherlands Nobel Biocare Netherlands Phone: +31 30 635 4949	New Zealand Nobel Biocare New Zealand Phone: +61 2 8064 5100 Cust. support: toll free 0800 441 657	Mexico Nobel Biocare Mexico Phone: +52 55 524 974 60	Ireland Nobel Biocare Ireland Phone: +353 1 800 677 306	Poland Nobel Biocare Poland Phone: +48 22 874 59 44 Cust. support: +48 22 874 59 45
---	---	---	---	--	--

Norway Nobel Biocare Norway Phone: +47 23 24 96 30	Poland Nobel Biocare Poland Phone: +48 22 874 59 44 Cust. support: +48 22 874 59 45	Portugal Nobel Biocare Portugal Phone: +351 22 374 73 50 Cust. support: toll free 800 300 100	Russia Nobel Biocare Russia Phone: +7 495 974 77 55	Singapore Nobel Biocare Singapore Phone: +65 6737 7967 Cust. support: toll free 800 53 93 28	South Africa Nobel Biocare South Africa Phone: +27 11 802 0112
---	--	--	--	---	---

Spain Nobel Biocare Spain Phone: +34 93 508 8800 Cust. support: toll free 900 850 008	Sweden Nobel Biocare Sweden Phone: +46 31 335 49 00 Cust. support: +46 31 335 49 10	Switzerland Nobel Biocare Switzerland Phone: +41 43 211 53 20	Taiwan Nobel Biocare Taiwan Phone: +886 2 27 93 99 33	United Kingdom Nobel Biocare UK Phone: +44 1895 430 650	USA Nobel Biocare USA Phone: +1 714 282 4800 Cust. support: +1 800 322 5501
--	--	--	--	--	--

ENGLISH

Important: Please read.

NobelProcera™ Abutment
Conical connection for NobelActive™
Internal connection for the following implant systems: NobelReplace™, Replace® Select and NobelSpeedy™ Replace.
External connection for the following implant systems: Brånemark System® and NobelSpeedy™ Groovy.
Depending on used NobelProcera™ software and scanner, abutments for other implant systems than from Nobel Biocare can be produced.

NobelProcera™ Abutments in both Zirconia and Titanium for: Straumann® Octagon 4.8mm, Camlog® Implant Systems 3.3, 3.8, 4.3, 5.0, 6.0mm.

NobelProcera™ Abutments in Titanium for: Astra Tech® Implant System 4.5ST, 5.0STmm. Astra Tech® Yellow, Aqua and Lilac. Straumann® Bone Level NC 3.3 and RC 4.1/4.8. Straumann® Octagon 6.5mm.

Disclaimer of liability: These prosthetic components are part of an overall concept and may only be used in conjunction with the associated original products according to the instructions and recommendation of Nobel Biocare. Use of products made by third parties in conjunction with Nobel Biocare prosthetic components will void any warranty or other obligation, express or implied, of Nobel Biocare. The user of Nobel Biocare prosthetic components has the duty to determine whether or not any product is suitable for the particular patient and circumstances. Nobel Biocare disclaims any liability, express or implied, and shall have no responsibility for any direct, indirect, punitive or other damages, arising out of or in connection with any errors in professional judgment or practice in the use of Nobel Biocare products. The clinician is also obliged to study the latest developments in regard to the Nobel Biocare products and their applications regularly. In cases of doubt, the clinician has to contact Nobel Biocare. Since the processing and surgical application of this product are under the control of the clinician, they are his/her responsibility. We assume no liability whatsoever for damage arising thereof.

Description: An individualized endosseous dental implant abutment directly connected to endosseous dental implant intended for use as an aid in prosthetic rehabilitation.

Indications: Indication range from single tooth to multiple teeth restorations using cement-retained supra-structures as well as screw-retained single tooth restorations.

Contraindication: NobelProcera™ Abutments Zirconia for NobelActive™ are not indicated for posterior use.

Instruction for clinician: We strongly recommend that clinicians, new as well as experienced implant users, always go through special training before undertaking a new treatment method. Nobel Biocare offers a wide range of courses for various levels of knowledge. For more info visit www.nobelbiocare.com <<http://www.nobelbiocare.com>>. Always work with an experienced colleague the first time you employ a new treatment method. Nobel Biocare has a global network of mentors for this purpose.

Procedural Precautions: Close cooperation between surgeon, restorative dentist and dental laboratory technician is essential for a successful implant treatment. It is especially important to achieve proper stress distribution through passive adaptation and fitting of the bridge to the implant abutments; adjusting occlusion to the opposing jaw; avoiding excessive transverse loading forces, particularly in immediate loading cases. If the prosthesis metal substructure is made of gold alloy, this should have a high gold content. Because of the small size of prosthetic components, care must be taken that they are not swallowed or aspirated by the patient.

Sterility: The NobelProcera™ Abutments co-packed with abutment screw are delivered non-sterile and should be cleaned and sterilized (if indicated) according to generally accepted procedures for dental laboratory work. Steam sterilize abutments and the prosthetic kit, if applicable for 5 min at 135°C/274°F.

Handling procedures

Clinical Procedure

- Place the impression coping onto the implant take an impression (**Illustration A**).
- Connect the healing abutment or temporary restoration.

Laboratory Procedure

- Produce a working model with removable gingival material.
- NobelProcera™ Abutment Wax-up:** Create the abutment in a pattern optical wax or resin, make sure this is made of a optical material – if not – surface needs to be coated with a conventional scanning spray (**B**). Use a temporary component (engaging) as a foundation for the abutment pattern. For NobelActive™, use a NobelProcera™ wax-up sleeve engaging NobelActive™ Internal. Scan the waxed abutment using the NobelProcera™ Scanner according to the tutorial found within the software.

Note: You must respect certain dimensions in the waxed abutment depending on which NobelProcera™ product you wish to fabricate:

- For ceramic abutments NP – maximum dimensions are 4mm x 12mm (Radius x Höhe).
- For ceramic abutments RP/WP – maximum dimensions are 5mm x 15mm (r x h).
- For titanium abutments NP – maximum dimensions are 5mm x 12mm (r x h).
- For titanium abutments RP/WP – maximum dimensions are 7,5mm x 15mm (r x h).

- NobelProcera™ Abutment CAD:** Insert an Abutment Position Locator into the replica (**C**). Scan the Abutment Position Locator using the NobelProcera™ Scanner and design the abutment using the NobelProcera™ software according to the tutorial found within the software.

After receiving the abutment from Nobel Biocare production facility.

- Finalizing procedures Zirconia Abutments:
 - If necessary, make minor adjustments with diamond bur or flex disc with fine grit size under low pressure and with copious water irrigation.
 - When designing/preparing NobelProcera™ Abutment Zirconia, make sure that the thickness of the ceramic material is at least 0.9mm. This thickness limit is applicable up to a height of 3mm above implant level (**D**). For NobelProcera™ Abutment NobelActive™, never modify the area of the abutment marked in dark. Do not modify the abutment below the dimensions shown (**E**).
 - Sandblast using max. one bar of pressure utilizing 110 µm aluminum oxide, at an approximate distance of 10 mm.
 - Clean in an ultrasonic unit.
- Finalizing procedures Titanium Abutments:
 - If necessary, make minor adjustments, with a carbide bur. Rubber wheels are used for the collar part.
 - Clean in an ultrasonic unit.
- Fabricate a crown or bridge with NobelProcera™ technique or with conventional casting technique. For single tooth screw-retained restorations it is also possible to veneer dental ceramic directly onto the abutment (**F**).

Clinical Procedure

- Connect the abutment onto the implant. It is recommended to verify the final abutment seating using radiographic imaging.
- Tighten the abutment according to table 1, using Screwdriver Unigrip™ and Manual Torque Wrench prosthetic (**G**).
Note: If final preparation is needed, make sure that the thickness of the ceramic material is at least 0.9mm and use the high-speed turbine with diamond drills and copious water irrigation. This thickness limit is applicable up to a height of 3mm above the implant level (**D**). For NobelProcera™ Abutment NobelActive™, never modify the area of the abutment marked in dark. Do not modify the abutment below the dimensions shown (**E**).
Note: Do not use temporary cement when cementing ceramic crowns and bridges due to increased risk of micro fractures.

NobelProcera™ Abutment Titanium and Zirconia screw chart (**table 2**).

For additional information on restorative and dental laboratory procedures please consult treatment guidelines available at www.nobelbiocare.com or request latest printed version from a Nobel Biocare representative.

Materials:
NobelProcera™ Abutments titanium - CP Titanium or Titanium alloy 90% Ti, 6% Al, 4% V.
NobelProcera™ Abutments zirconia - Yttria-stabilized Zirconiumoxide.
Abutment screws - Titanium alloy 90% Ti, 6% Al, 4% V.

Manufacturer: Nobel Biocare AB, Box 5190, SE-402 26 Göteborg, Sweden.
Phone: +46 31 81 88 00. Fax: +46 31 16 31 52. www.nobelbiocare.com

Caution: The caution text "Federal (USA) law restricts the sale of this device to, or on the order of, a licensed physician or dentist" is shown on labels with "Rx Only".

Straumann® is a registered trademark of Institut Straumann AG.
Camlog® is a trademark of Camlog Biotechnologies Group.
Astra Tech® is a trademark of Astra Tech Group.

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