NobelProcera™ Crown and Bridge

Crown and Bridge in Alumina.

73291 West 1105 Printed in Sweden. Date of issue 2011-03-29.

All rights reserved. Nobel Biocare, the Nobel Biocare logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Nobel Biocare. Cerabien is a registered trademark of Noritake Co., Limited. VITA VM7 is a registered trademark of VITA Zahnfabrik H. Rauter GmbH & Co.KG. Disclaimer: Some products may not be regulatory cleared/released for sale in all markets. Please contact the local Nobel Biocare sales office for current product assortment and availability.

ENGLISH

Important: Please read.

Disclaimer of liability: This prosthetic component is part of an overall concept and should only be used in conjunction with the recommended products according to the instructions and recommendation of Nobel Biocare. Use of non-recommended products made by third parties in conjunction with Nobel Biocare prosthetic components will void any warranty or other obligation, express or implied, of Nobel Biocare. The user of Nobel Biocare prosthetic components has the duty to determine whether or not any product is suitable for the particular patient and circumstances. Nobel Biocare disclaims any liability, express or implied, and shall have no responsibility for any direct, indirect, punitive or other damages, arising out of or in connection with any errors in professional judgment or practice in the use of Nobel Biocare products. The clinician is also obliged to study the latest developments in regard to the Nobel Biocare products and their applications regularly. In cases of doubt, the clinician has to contact Nobel Biocare. Since the processing and application of this product are under the control of the clinician, they are his/her responsibility. We assume no liability whatsoever for damage arising thereof. Failure to observe the indicated limitations of use and working steps may result in failure. Please note that not all products may have been licenced in accordance with Canadian law.

Description: An individualized coping or bridge cement retained solution on natural teeth or abutments, for use as the core structure for a crown and bridge.

Materials: NobelProcera™ Alumina – Biocompatible densely sintered high purity α-alumina (Al₂O₃).

Indications:

NobelProcera™ Crown Alumina:

- Cement retained solution on teeth and abutments.
- Applicable for all indications, however best used in esthetic region due to its translucent properties.
- NobelProcera™ Crown Alumina is available in the thicknesses:
 - 0.4 mm, is indicated for use in the esthetic region, premolar to premolar.
 - 0.6 mm for all indications.

Contraindication:

NobelProcera™ Crown Alumina:

- Bruxism.

Indications:

NobelProcera™ Bridge Alumina:

- Cement retained solution on teeth and abutments.
- Esthetic region, premolar to premolar.
- Bridge of 2–4 units.
- Connector requirements; minimum cross-sectional area of 6.0mm² and minimum height 3.0mm and minimum width 2.5mm.

Contraindication:

NobelProcera™ Bridge Alumina:

- No cantilevers and a maximum of one pontic between supporting teeth/abutments.
- Bruxism.

Procedural Precautions: Close cooperation between restorative dentist and dental laboratory technician is essential for a successful treatment. It is especially important to achieve proper stress distribution through adaptation and fitting of the crown by adjusting occlusion to the opposing jaw. In addition avoid excessive transverse loading forces, particularly in immediate loading cases. Because of the small size of prosthetic components, care must be taken that they are not swallowed or aspirated by the patient. If the restoration requires modification, pay attention to the marginal area and occlusal surfaces; a minimum thickness of 0.4 mm should be kept for copings and 0.6 mm for bridges. If modifying the restoration, use appropriate protection equipment to avoid inhale the NobelProcera™ Alumina dust. If a try-in on the patient is carried out, we recommend cleaning the restoration under running water or with alcohol.

Sterility: The coping or bridge is delivered non-sterile and should be disinfected with a dedicated dis-infection solution (e.g 70% ethanol) before try-in and again after try-in directly before cementation.

Handling procedures

Clinical Procedure

- Preparation of tooth:**
 - Eliminate sharp edges, undercuts, and grooves.
 - Provide sufficient reduction (1,2 mm to 1,5 mm).
 - Provide adequate space for the crown.
 - Avoid sharp angles on the occlusal surface.
 - Avoid creating a deep fossa/cavity.
 - Provide sufficient (1,5 mm to 2 mm) occlusal/incisal reduction.
 - Avoid preparations that are excessively tapered or too close to parallel. The ideal total occlusal convergence is 6–10°.
- Impression:**
 - Make an ordinary crown and bridge impression using your preferred impression material and method.
 - Send the impression to the laboratory.*Procedure after the restoration has been delivered from the dental laboratory*
- Cementation:** Gently seat the restoration on the tooth and check both the occlusion and the interproximal contacts. The restoration should be in light occlusion. Excursive contact should be minimal.
- Removal:** To avoid overheating cross-cutting burs with copious irrigation during preparation should be used, if rotating instruments are used for removal.

Laboratory Procedure

- Production:** A model is made from the impression and scanned for producing a NobelProcera™ Crown or bridge with the NobelProcera™ Scanner and the NobelProcera™ software according to the tutorial found within the software. The framework should be designed with optimum veneering material support. The coping or bridge is milled in a NobelProcera™ production site.
- Veneering:**
 - All commercially available veneering materials designed for use with Alumina with a CTE of 7 are suitable. A veneering option is VITA VM® 7 or Noritake Cerabien.
 - For long-term clinical success follow the recommendations and handling instructions of the veneering material manufacturer.

--	--	--	--	--	--

Manufacturer: Nobel Biocare AB, Box 5190, SE-402 26 Göteborg, Sweden. Phone: +46 31 81 88 00. Fax: +46 31 16 31 52. www.nobelbiocare.com

Caution: The caution text "Federal (USA) law restricts the sale of this device to, or on the order of, a licensed physician or dentist" is shown on labels with "Rx Only".

	Do not re-use. Re-use could result in local infection or transmission of infectious diseases.		
--	---	--	--

Argentina Nobel Biocare Argentina Phone: +54 11 4825 9696 Cust. support: toll free 0800 800 66235	Belgium Nobel Biocare Belgium Phone: +32 2 467 41 70	China Nobel Biocare China Phone: +86 21 5206 6655 Cust. support: +86 21 5206 0974	France Nobel Biocare France Phone: +33 1 49 20 00 30	Germany Nobel Biocare Germany Phone: +49 221 500 585 590	Hong Kong Nobel Biocare Hong Kong Cust. support: +852 2823 8926	India Nobel Biocare India Phone: +91 22 6751 9999 Cust. support: toll free 1 800 22 9998	Indonesia Nobel Biocare Indonesia Phone: +62 21 5727 7777 Cust. support: toll free 1 800 804 597	Ireland Nobel Biocare Ireland Phone: toll free 1800 677 306	Italy Nobel Biocare Italy Phone: +39 039 683 61 Cust. support: toll free 800 53 93 28	Japan Nobel Biocare Japan Phone: +81 3 6717 6191	Lithuania Nobel Biocare Lithuania Phone: +370 5 268 3448 Cust. support: toll free 880 01 23 24	Mexico Nobel Biocare Mexico Phone: +52 55 244 974 60	Middle East Nobel Biocare Middle East Phone: +48 22 874 5951	Netherlands Nobel Biocare Netherlands Phone: +31 30 635 4949	New Zealand Nobel Biocare New Zealand Phone: +61 2 8064 5100 Cust. support: toll free 0800 441 657	Norway Nobel Biocare Norway Phone: +47 23 24 98 30	Poland Nobel Biocare Poland Phone: +48 22 874 59 44 Cust. support: +48 22 874 59 45	Portugal Nobel Biocare Portugal Phone: +351 22 374 73 50 Cust. support: toll free 800 300 100	Russia Nobel Biocare Russia Phone: +7 495 974 77 55 Cust. support: toll free 8 800 250 77 55	Singapore Nobel Biocare Singapore Phone: +65 6737 7967 Cust. support: +65 6737 7967	South Africa Nobel Biocare South Africa Phone: +27 11 802 0112	Spain Nobel Biocare Spain Phone: +34 93 508 8800 Cust. support: toll free 900 850 008	Sweden Nobel Biocare Sweden Phone: +46 31 335 49 00 Cust. support: +46 31 335 49 10	Switzerland Nobel Biocare Switzerland Phone: +41 43 211 53 20	Taiwan Nobel Biocare Taiwan Phone: +886 2 2793 9933	United Kingdom Nobel Biocare UK Phone: +44 1895 430 650	USA Nobel Biocare USA Phone: +1 714 282 4800 Cust. support: +1 800 322 5001
--	---	--	---	---	---	---	---	--	--	---	---	---	---	---	---	---	--	--	---	--	---	--	--	--	--	--	--

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Fabricatore: Nobel Biocare AB, Box 5190, SE-402 26 Göteborg, Svezia