
NobelProcera™ Crown and Bridge

Crown and Bridge in Zirconia.

INSTRUCTIONS FOR USE

72835 ML West 1103 Printed in Sweden. Date of issue 2011-02-09.
All rights reserved. Nobel Biocare, the Nobel Biocare logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Nobel Biocare.
VITA VM9 and VITA PM9 are registered trademarks of VITA Zahnfabrik H. Rauter GmbH & Co.KG. IPS e.max Ceram and IPS e.max ZirPress are registered trademarks of Ivoclar Vivadent AG. Product images in this folder are not necessarily to scale.

Contraindications:
– Bruxism
– Cases where the mesial/distal cantilever have a length of more than one unit.
– External radius less than 0.6mm.

Procedural Precautions: Close cooperation between restorative dentist and dental laboratory technician is essential for a successful treatment. It is especially important to achieve proper stress distribution through adaptation and fitting of the crown by adjusting occlusion to the opposing jaw. In addition avoid excessive transverse loading forces, particularly in immediate loading cases. Because of the small size of prosthetic components, care must be taken that they are not swallowed or aspirated by the patient. If the restoration requires modification, pay attention to the marginal area and occlusal surfaces; a minimum thickness of 0.4 mm should be kept. If modifying the restoration, use appropriate protection equipment to avoid inhale the NobelProcera™ Zirconia dust. If a try-in on the patient is carried out, we recommend cleaning the restoration under running water or with alcohol.

Sterility: The coping or bridge is delivered non-sterile and should be disinfected with a dedicated disinfection solution (e.g. 70% ethanol) before try-in and again after try-in directly before cementation.

Handling procedures

Clinical Procedure

- Preparation of tooth:**
 - Eliminate sharp edges, undercuts, and grooves.
 - Provide sufficient reduction (1,2 mm to 1,5 mm).
 - Provide adequate space for the crown.
 - Avoid sharp angles on the occlusal surface.
 - Avoid creating a deep fossa/cavity.
 - Provide sufficient (1,5 mm to 2 mm) occlusal/incisal reduction.
 - Avoid preparations that are excessively tapered or too close to parallel.
The ideal total occlusal convergence is 6°–10°.
- Impression:**
 - Make an ordinary crown and bridge impression using your preferred impression material and method.
 - Send the impression to the laboratory.
Procedure after the restoration has been delivered from the dental laboratory
- Cementation:** Gently seat the restoration on the tooth and check both the occlusion and the interproximal contacts. The restoration should be in light occlusion. Excursive contact should be minimal.
- Removal:** To avoid overheating cross-cutting burs with copious irrigation during preparation should be used, if rotating instruments are used for removal.

Laboratory Procedure

- Production:** A model is made from the impression and scanned for producing a NobelProcera™ Crown or bridge with the NobelProcera™ Scanner and the NobelProcera™ software according to the tutorial found within the software. The framework should be designed with optimum veneering material support. The coping or bridge is milled in a NobelProcera™ production site.
- Veneering:** After production, the crown is veneered with veneering material designed for use with zirconium with a CTE of 10.4. A veneering option is Ivoclar Vivadent IPS e.max Ceram and e.max ZirPress or VITA VM9 and PM9.

Manufacturer: Nobel Biocare AB, Box 5190, SE-402 26 Göteborg, Sweden.
Phone: +46 31 81 88 00. Fax: +46 31 16 31 52. www.nobelbiocare.com

Caution: The caution text "Federal (USA) law restricts the sale of this device to, or on the order of, a licensed physician or dentist" is shown on labels with "Rx Only".

Do not re-use. Re-use could result in local infection or transmission of infectious diseases.

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ENGLISH

Important: Please read.

Disclaimer of liability: This prosthetic component is part of an overall concept and should only be used in conjunction with the recommended products according to the instructions and recommendation of Nobel Biocare. Use of non-recommended products made by third parties in conjunction with Nobel Biocare prosthetic components will void any warranty or other obligation, express or implied, of Nobel Biocare. The user of Nobel Biocare prosthetic components has the duty to determine whether or not any product is suitable for the particular patient and circumstances. Nobel Biocare disclaims any liability, express or implied, and shall have no responsibility for any direct, indirect, punitive or other damages, arising out of or in connection with any errors in professional judgment or practice in the use of Nobel Biocare products. The clinician is also obliged to study the latest developments in regard to the Nobel Biocare products and their applications regularly. In cases of doubt, the clinician has to contact Nobel Biocare. Since the processing and application of this product are under the control of the clinician, they are his/her responsibility. We assume no liability whatsoever for damage arising therefrom. Failure to observe the indicated limitations of use and working steps may result in failure. Please note that not all products may have been licensed in accordance with Canadian law.

Description: An individualized coping or bridge cement retained solution on natural teeth or abutments, for use as the core structure for a crown and bridge.


Materials: NobelProcera™ Zirconia – Biocompatible yttria-stabilized tetragonal zirconia (Y-TZP).

Indications:

- Anterior or posterior coping or bridge for functional and esthetic replacement of teeth.
- Tooth- and implant-supported restorations.
- NobelProcera™ Crown Zirconia, 0.4mm, is indicated for use in the esthetic region, premolar to premolar.
- 0.7 mm for all locations in the mouth.

- For a bridge span of 2–14 units, designed to fit into a disc, 60 mm in diameter and 20 mm in height.
- 0.6 mm for NobelProcera™ Zirconia bridges, any position in the mouth.
- The connector dimension of a multi-unit framework depends on the distance between the supporting teeth (see table below). The design is supported by real-time feedback through warning function in the NobelProcera™ Software.

Type – any position	Maximum length	Minimum connector and cross section
	a=Arc length [mm] e=Extension length [mm]	height x width [mm] / area [mm ²]
Free hanging arc	a ≤ 21,0	3,0 x 2,5 / 6,0
Free hanging arc	21,0 < a ≤ 35,0	4,0 x 3,0 / 9,4
Extension	e ≤ 10,0	4,0 x 3,0 / 9,4

Contraindications:
– Bruxism
– Cases where the mesial/distal cantilever have a length of more than one unit.
– External radius less than 0.6mm.

Procedural Precautions: Close cooperation between restorative dentist and dental laboratory technician is essential for a successful treatment. It is especially important to achieve proper stress distribution through adaptation and fitting of the crown by adjusting occlusion to the opposing jaw. In addition avoid excessive transverse loading forces, particularly in immediate loading cases. Because of the small size of prosthetic components, care must be taken that they are not swallowed or aspirated by the patient. If the restoration requires modification, pay attention to the marginal area and occlusal surfaces; a minimum thickness of 0.4 mm should be kept. If modifying the restoration, use appropriate protection equipment to avoid inhale the NobelProcera™ Zirconia dust. If a try-in on the patient is carried out, we recommend cleaning the restoration under running water or with alcohol.

Sterility: The coping or bridge is delivered non-sterile and should be disinfected with a dedicated disinfection solution (e.g. 70% ethanol) before try-in and again after try-in directly before cementation.

Handling procedures

Clinical Procedure

- Preparation of tooth:**
 - Eliminate sharp edges, undercuts, and grooves.
 - Provide sufficient reduction (1,2 mm to 1,5 mm).
 - Provide adequate space for the crown.
 - Avoid sharp angles on the occlusal surface.
 - Avoid creating a deep fossa/cavity.
 - Provide sufficient (1,5 mm to 2 mm) occlusal/incisal reduction.
 - Avoid preparations that are excessively tapered or too close to parallel.
The ideal total occlusal convergence is 6°–10°.
- Impression:**
 - Make an ordinary crown and bridge impression using your preferred impression material and method.
 - Send the impression to the laboratory.
Procedure after the restoration has been delivered from the dental laboratory
- Cementation:** Gently seat the restoration on the tooth and check both the occlusion and the interproximal contacts. The restoration should be in light occlusion. Excursive contact should be minimal.
- Removal:** To avoid overheating cross-cutting burs with copious irrigation during preparation should be used, if rotating instruments are used for removal.

Laboratory Procedure

- Production:** A model is made from the impression and scanned for producing a NobelProcera™ Crown or bridge with the NobelProcera™ Scanner and the NobelProcera™ software according to the tutorial found within the software. The framework should be designed with optimum veneering material support. The coping or bridge is milled in a NobelProcera™ production site.
- Veneering:** After production, the crown is veneered with veneering material designed for use with zirconium with a CTE of 10.4. A veneering option is Ivoclar Vivadent IPS e.max Ceram and e.max ZirPress or VITA VM9 and PM9.

Manufacturer: Nobel Biocare AB, Box 5190, SE-402 26 Göteborg, Sweden.
Phone: +46 31 81 88 00. Fax: +46 31 16 31 52. www.nobelbiocare.com

Caution: The caution text "Federal (USA) law restricts the sale of this device to, or on the order of, a licensed physician or dentist" is shown on labels with "Rx Only".

Do not re-use. Re-use could result in local infection or transmission of infectious diseases.

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