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English
Figure: Straumann® CAD/CAM abutment RC, ceramic (to exemplify the entire Straumann® CAD/CAM ceramic abutment portfolio)

Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed dental professional. [Rx Only]

Product description

Straumann® CAD/CAM abutments are used for the restoration of Straumann® dental implants of different types, endosteal diameters, lengths and platforms. The Straumann® CAD/CAM abutments allow for individual customization regarding function and esthetics. Straumann® CAD/CAM abutments are designed by the customer either by means of a traditional waxup abutment that is subsequently scanned or by scanning of the intraoral situation and designing of the shape by the CAD module of the integrated software. The design data is then sent to Straumann where the fabrication of the customized abutment is carried out. These instructions for use are valid for Straumann® CAD/CAM ceramic abutments of the following prosthetic lines: synOcta® (RN), CrossFit™ (NC and RC) and Narrow Neck (NN).

Indications for use

Straumann® CAD/CAM ceramic abutments are intended to be placed into Straumann dental implants to provide support for prosthetic reconstructions such as crowns, bridges and overdentures.

The **Straumann® CAD/CAM ceramic abutments** are indicated for single-tooth replacement and multiple tooth restorations. The prosthetic restoration can be cemented or directly veneered/screw-retained.

Contraindications

Allergies to materials used, which may include any or all of the following: zirconium dioxide (ZrO₂), titanium (Ti), titanium alloy (Ti-6Al-7Nb or TAN).

Important

Straumann® CAD/CAM ceramic abutments supplied must not be ground or polished after delivery. The handling and use of the product are not within the control of Straumann, and are the responsibility of the user. All liability for loss or damage is excluded.

General warnings & precautions

- Our products must be secured to prevent aspiration during intraoral use or insertion into the implant
- Aseptic handling after sterilization procedure is essential

Instructions provided are insufficient to serve as the only means for placing the Straumann abutments and related components. Only those clinicians thoroughly trained in dental implantology should be placing these devices. Placing abutments and related components without proper training may lead to failure of the component or abutment. Abutment failure may lead to implant removal or other complications.

For further detail, see additional Straumann brochures identified below.

Failure to follow the procedures outlined in these instructions may harm the patient and/or lead to any or all of the following complications:

- Aspiration of components;
- Damage to the implant, abutment, components or tooling;
- Loosening of the abutment or other components;
- Improper final restoration or malfunction of the crown, bridge, or other final prosthesis;
- Impairment of the patient's chewing function;
- Failure of the implant; and/or
- Removal of the implant.

Provisional restorations must always be placed out of occlusion.

Provisional cement, cement or any other material used for attaching prosthetic components to others must be processed as specified by the manufacturer.

Compatibility information

Straumann implants and the prosthetic lines synOcta®, CrossFit™ and Narrow Neck are available in a variety of configurations to meet your clinical needs. The label on each product uses abbreviations to help you identify whether a particular abutment or coping is compatible with the implant that you are restoring.

Straumann's Bone Level (BL) implant line uses the terms RC (Regular CrossFit™) for 4.1 and 4.8 mm implant diameters or NC (Narrow CrossFit™) for 3.3 mm implant diameters. RC BL implants are only compatible with components marked RC. NC BL implants are only compatible with components marked NC.

Straumann's Standard (S), Standard Plus (SP) and Tapered Effect (TE) implant lines use the term RN (Regular Neck) for implants with 4.8 mm restorative platforms. RN S, RN SP, and RN TE implants are only compatible with components marked RN.

The term NN (Narrow Neck) describes Straumann's 3.5 mm restorative platform of the Standard Plus Implant Narrow Neck. SP NN implants are only compatible with components marked NN.

Make sure that you use only original Straumann parts with the corresponding connection for restoring of a Straumann implant.

The following table summarizes the information above:

Implant type	Connection type	Compatible parts
BL (Bone Level)	RC (Regular CrossFit™)	Parts labelled RC
	NC (Narrow CrossFit™)	Parts labelled NC
S (Standard)	RN (Regular Neck)	Parts labelled RN
SP (Standard Plus)	TE (Tapered Effect)	
NN (Standard Plus/ Narrow Neck)	NN (Narrow Neck)	Parts labelled NN

Note: You must only restore an implant with abutments and copings that are compatible with that implant. Failure to follow this instruction may harm the patient and/or result in damage to the implant, components or tooling.

Use and handling of Straumann prosthetic ceramic components

1. Sterilization

Straumann® abutments and components are not sterile when delivered. Straumann recommends the following procedure for sterilization prior to use:

Component	Material	Method	Conditions
Abutment	ZrO ₂	Dry heat	160 °C (320 °F), 4 h
Screw, RN synOcta® 1.5 abutment	Ti, Ti alloy	Autoclave, moist heat	134 °C (273 °F), 18 min

Note: Parts that have been modified or altered in any way may require a different sterilisation procedure.
Note: ZrO₂ abutments must neither be autoclaved nor chemiclaved.

Caution: Straumann® CAD/CAM ceramic abutments must not be ground or polished after delivery. Grounding or polishing Straumann® CAD/CAM ceramic abutments may cause the product to fail.

2. Straumann secondary parts are inserted into the implant without applying cement. Straumann® CAD/CAM abutments are inserted and tightened using a Straumann® SCS screwdriver.

3. Components that are secured to the implant with a premounted screw can be brought to the mouth using the SCS screwdriver. All other components should be handled fully using tweezers.

4. Insert the abutment or the component into the connection.

5. Make sure that the retentive elements of the implant abutment connection are properly aligned for all secondary parts such as abutments as described in the table below.

Implant type	Implant/Secondary part alignment
BL (Bone Level)	Internal grooves
S (Standard)	Internal octagon
SP (Standard Plus)	
TE (Tapered Effect)	
NN (Narrow Neck)	External octagon

6. Secure the mounted secondary part to the implant with the appropriate screw using the Straumann SCS screwdriver, ratchet and torque control device. Secure the mounted **RN Straumann® CAD/CAM ceramic abutment** to the RN synOcta 1.5 screw-retained abutment using the Straumann SCS screwdriver, ratchet and torque control device.

Straumann secondary parts are inserted into the implant without applying cement. The **RN Straumann® CAD/CAM ceramic abutment** is mounted on the RN synOcta 1.5 screw-retained abutment without applying cement.

7. Use the ratchet to torque the screw until the appropriate tightening torque is reached [refer to the table below]. For stabilization, use the Straumann holding key. Ensure the arrow of the torque control device is pointing clockwise [towards the torque bar with the leadrop]. If it's not, simply pull the arrow out, flip over and push it back in. Using one hand on the holding key, use one finger from the other hand to push the leadrop of the torque control device. Apply pressure to the torque bar until the required torque is reached. After reaching the indicated tightening torque, return the torque bar to its starting position. Lift and remove the holding key, ratchet with torque control device, and the driver.

Device type	Tightening torque
Abutments (permanent)	35 Ncm
RN CAD/CAM abutment on RN synOcta® 1.5 abutment (permanent)	35 Ncm
Components on implant analogs	Hand tight

Warning

Torques greater than 35 Ncm may result in failure of the abutment and/or implant. Torque values less than the recommended values may result in loosening of the abutment, which may lead to abutment and/or implant failure.

Cemented restorations and screw-retained directly veneered mesostructures

A: Cementing the crown or bridge

Before cementing the crown or bridge, the screw opening(s) is [are] sealed off with wax or gutta-percha. This makes it possible to release the screw again if required. Then final cementing of the crown to the mesostructure can take place.

B: Closure of the screw channel with directly veneered mesostructure
To enable the occlusal screw to be released again, if required, at a later point in time, the screw head is covered with wax or gutta-percha and the screw channel is sealed off with a temporary veneering material (e.g. composite).

Firing parameters for direct ceramic veneering (for B)

The coefficient of thermal expansion of the abutment is in the 10.5 x 10⁻⁶/K (25–500 °C, 77–932 °F) range. According to this guide, the product can be veneered with suitable dental ceramics available on the market. For more detailed information about direct ceramic veneering, follow the instructions for use provided by the manufacturer of the veneering porcelain. **Mechanical finishing of the abutment (e.g. sandblasting) is not permitted.**

General instructions

To ensure correct transfer of the position of the Straumann® CAD/CAM ceramic abutment from the master cast to the patient, an individual index can be fabricated on the cast. In the case of single crowns, the index is secured with support from the adjacent teeth, and in the case of bridges the abutments are splinted to one another.

Important

The screw opening must not be covered with plastic. Ensure that no plastic gets into the interior of the abutment.

Further information

For additional information about the use of Straumann products, call Straumann's customer service department. The following brochures provide additional guidance on the use of the Straumann® prosthetic components and instruments.

- “Straumann® CAD/CAM abutments – Basic information on the dental lab procedures“, Art. No. 152.099
- “Crown and bridge restorations – Straumann® synOcta® Prosthetic System“, Art. No. 152.255
- “Basic information on the prosthetic procedures – Straumann® Bone Level Implant line“, Art. No. 152.810

Disclaimer of liability

The Straumann® dental implants and other Straumann products are part of an overall concept and may be used only in conjunction with the corresponding original components and instruments according to Institut Straumann AG's instructions and recommendations.

Use of products made by third parties, which are not distributed through Institut Straumann AG, in conjunction with the Straumann Dental Implant System will void any warranty or other obligation, express or implied, of Straumann.

Instructions as to the application of our products take place verbally, in writing, by electronic media or in hands-on trainings corresponding to the state of the art at the time of introduction of the product.

The user of Straumann products is responsible for determining whether or not any product is suitable for a

 straumann
Instruções de uso: Pilares Straumann® monochoni CADCAM, cerâmico
Bruksanvisning: Straumann® CADCAM distanser, keramik
Gebrauchsanweisung: Straumann® CADCAM-abutments, keramik
Brugsanvisning: Straumann® CADCAM abutments, keramik
Инструкция по применению: опорные зубы Straumann® CADCAM, керамика
Fabricante/Fabrikante/Tillverkare/Producent/Изготовитель: Institut Straumann AG, CH-4002 Basel/Switzerland, www.straumann.com

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Copyright und Warenzeichen