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English

Figure: Straumann® CAD/CAM abutment RC, titanium (to exemplify the entire Strau-mann® CAD/CAM titanium abutment portfolio)

Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed dental professional. [Rx Only]

Product description

Straumann® CAD/CAM abutments are used for the restoration of Straumann® dental implants of different types, endosteal diameters, lengths and platforms. The Straumann® CAD/CAM abutments allow for individual customization regarding function and esthetics. Straumann® CAD/CAM abutments are designed by the customer either by means of a traditional waxup abutment that is subsequently scanned or by scanning of the intraoral situation and designing of the shape by the CAD module of the integrated software. The design data is then sent to Straumann where the fabrication of the customized abutment is carried out. These instructions for use are valid for Straumann® CAD/CAM titanium abutments of the following prosthetic lines: synOcta® (RN and WN), CrossFit™ (NC and RC) and Narrow Neck (NN).

Indications for use

Straumann® CAD/CAM titanium abutments are intended to be placed into Straumann dental implants to provide support for prosthetic reconstructions such as crowns, bridges and overdentures.

The **Straumann® CAD/CAM titanium abutments** are indicated for cemented single-tooth replacements and multiple tooth restorations.

Contraindications

Allergies to materials used, which may include any or all of the following: titanium (Ti), titanium alloy (Ti-6Al-7Nb or TAN).

Important

Straumann® CAD/CAM titanium abutments must not be directly veneered with ceramic.

Straumann® CAD/CAM titanium abutments supplied must not be ground or polished after delivery. The handling and use of the product are not within the control of Straumann, and are the responsibility of the user. All liability for loss or damage is excluded.

General warnings & precautions

- Our products must be secured to prevent aspiration during intraoral use or insertion into the implant.
- Titanium secondary parts must not be directly veneered with ceramic. Failure to follow this instruction could result in damage to the abutment (heating up to 300 °C /572 °F) or more decreases the strength of titanium parts).
- Aseptic handling after sterilization procedure is essential

Instructions provided are insufficient to serve as the only means for placing the Strau-mann abutments and related components. Only those clinicians thoroughly trained in dental implantology should be placing these devices. Placing abutments and related components without proper training may lead to failure of the component or abutment. Abutment failure may lead to implant removal or other complications.

For further details, see additional Straumann brochures identified below.

Failure to follow the procedures outlined in these instructions may harm the patient and/ or lead to any or all of the following complications:

- Aspiration of components;
- Damage to the implant, abutment or other components;
- Loosening of the abutment or other components;
- Improper final restoration or malfunction of the crown, bridge, or other final prosthetic;
- Impairment of the patient's chewing function;
- Failure of the implant; and/or
- Removal of the implant.

Provisional restorations must always be placed out of occlusion.

Provisional cement, cement or any other material used for attaching prosthetic components to others must be processed as specified by the manufacturer.

Compatibility information

Straumann implants and the prosthetic lines synOcta®, CrossFit™ and Narrow Neck are available in a variety of configurations to meet your clinical needs. The label on each product uses abbreviations to help you identify whether a particular abutment or coping is compatible with the implant that you are restoring.

Straumann's Bone Level (BL) implant line uses the terms RC (Regular CrossFit™) for 4.1 and 4.8 mm implant diameters or NC (Narrow CrossFit™) for 3.3 mm implant diam-eters. RC BL implants are only compatible with components marked RC. NC BL implants are only compatible with components marked NC.

Straumann's Standard (S), Standard Plus (SP) and Tapered Effect (TE) implant lines use the terms RN (Regular Neck) for implants with 4.8 mm restorative platforms and WN (Wide Neck) for implants with 6.5 mm restorative platforms. RN S, RN SP, and RN TE implants are only compatible with components marked RN. WN S, WVN SP and WVN TE implants are only compatible with components marked WN.

The term NN (Narrow Neck) describes Straumann's 3.5 mm restorative platform of the Standard Plus Implant Narrow Neck. SP NN implants are only compatible with components marked NN.

Make sure that you use only original Straumann parts with the corresponding connection for restoring of a Straumann implant.

The following table summarizes the information above:

Implant type	Connection type	Compatible parts
BL (Bone Level)	RC (Regular CrossFit™)	Parts labelled RC
	NC (Narrow CrossFit™)	Parts labelled NC
S (Standard), SP (Standard Plus), TE (Tapered Effect)	WN (Wide Neck)	Parts labelled WN
	RN (Regular Neck)	Parts labelled RN
NN (Standard Plus Narrow Neck)	NN (Narrow Neck)	Parts labelled NN

Note: You must only restore an implant with abutments and copings that are compatible with that implant. Failure to follow this instruction may harm the patient and/ or result in damage to the implant, components or tooling.

Use and handling of Straumann prosthetic titanium components

1. Sterilization

Straumann® abutments and components are not sterile when delivered. Straumann rec-ommends the following procedure for sterilization prior to use:

Material	Method	Conditions
Ti, Ti alloy	Autoclave, moist heat	134 °C [273 °F], 18 min

Note: Parts that have been modified or altered in any way may require a diffe-

rent sterilisation procedure.

Caution: Straumann® CAD/CAM titanium abutments must not be ground or polished after delivery. Grinding or polishing Straumann® CAD/CAM titanium abutments may cause the product to fail.

2. Straumann secondary parts are inserted into the implant without applying cement. Straumann® CAD/CAM abutments are inserted and tightened using a Straumann® SCS screwdriver.

3. Make sure that the retentive elements of the implant abutment connection are prop-erly aligned for all secondary parts such as abutments as described in the table below.

Implant type	Implant/Secondary part alignment
BL (Bone Level)	Internal grooves
	Internal octagon
S (Standard), SP (Standard Plus), TE (Tapered Effect)	Internal octagon
	External octagon

4. Secure the mounted secondary part to the implant with the appropriate screw using the Straumann SCS screwdriver, ratchet and torque control device.

5. Use the ratchet to torque the screw until the appropriate tightening torque is reached (refer to the table below). For stabilization, use the Straumann holding key. Ensure the arrow of the torque control device is pointing clockwise (towards the torque bar with the leadrop). If it's not, simply pull the arrow out, flip over and push it back in. Using one hand on the holding key, use one finger from the other hand to push the leadrop only of the torque control device. Apply pressure to the torque bar until the required torque is reached. After reaching the indicated tightening torque, return the torque bar to its starting position. Lift and remove the holding key, ratchet with torque control device, and the driver.

Device type	Tightening torque
Abutments (permanent)	35 Ncm
Components on implant analogs	Hand tight

Warning

Torques greater than 35 Ncm may result in failure of the abutment and/or implant. Torque values less than the recommended values may result in loosening of the abutment, which may lead to abutment and/ or implant failure.

Cemented crown or bridge

Before cementing the crown or bridge, the screw opening(s) is (are) sealed off with wax or gutta-percha. This makes it possible to release the screw again if required. Then final cementing of the crown to the mesostructure can take place.

General instructions

To ensure correct transfer of the position of the Straumann® CAD/CAM titanium abutment from the master cast to the patient, an individual index can be fabricated on the cast. In the case of single crowns, the index is secured with support from the adjacent teeth, and in the case of bridges the abutments are splinted to one another.

Important:

The screw opening must not be covered with plastic. Ensure that no plastic gets into the interior of the abutment.

Further information

For additional information about the use of Straumann products, call Straumann's cus-tomer service department. The following brochures provide additional guidance on the use of the Straumann® prosthetic components and instruments.

- “Straumann® CAD/CAM abutments – Basic information on the dental lab proce-dures”, Art. No. 152.099
- “Crown and bridge restorations – Straumann® synOcta® Prosthetic system”, Art. No.152.255
- “Basic information on the prosthetic procedures – Straumann® Bone Level Implant line”, Art. No. 152.810

Disclaimer of liability

The Straumann® dental implants and other Straumann products are part of an overall concept and may be used only in conjunction with the corresponding original com-ponents and instruments according to Institut Straumann AG's instructions and recom-mendations.

Use of products made by third parties, which are not distributed through Institut Strau-mann AG, in conjunction with the Straumann Dental Implant System will void any warranty or other obligation, express or implied, of Straumann.

Instructions as to the application of our products take place verbally, in writing, by electronic media or in hands-on trainings corresponding to the state of the art at the time of introduction of the product.

The user of Straumann products is responsible for determining whether or not any product is suitable for a particular patient and circumstances. Straumann disclaims any liability, express or implied, and bears no responsibility for any direct, indirect, punitive or other damages, arising out of or in connection with any errors in professional judg-ment or practice in the use or installation of Straumann products.

The user is also obliged to study the latest developments of the Straumann Dental Implant System and their applications regularly.

Please note

The descriptions contained in this document are not sufficient for immediate use of the Straumann Dental Implant System. Knowledge of dental implantology and instruction in the handling of the Straumann Dental Implant System provided by an operator with the relevant experience are always necessary.

Supply, availability

Some items in the Straumann® Dental Implant System are not available in all countries.

Copyright and trademarks

Straumann® and/or other products and logos from Straumann® that are mentioned here are trademarks or registered trademarks of Straumann Holding AG and/or its affiliates.

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