72219 C ML West 1163 Printed in Sweden. Date of issue 2010-09-02. All rights reserved. Nobel Biocare, the Nobel Biocare logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Nobel Biocare. IPS e.max® is a trademark of Ivoclar Vivadent AG. Product images in this folder are not necessarily to scale.

ENGLISH

Important: Please read.

Disclaimer of liability: These prosthetic components are part of an overall concept and may only be used according to the instructions and recommendation of Nobel Biocare. Non-recommended use of products made by third parties in conjunction with Nobel Biocare prosthetic components will void any warranty or other obligation, express or implied, of Nobel Biocare. The user of Nobel Biocare prosthetic components has the duty to determine whether or not any product is suitable for the particular patient, and circumstances. Nobel Biocare disclaims any liability, express or implied, and shall have no responsibility for any direct, indirect, punitive or other damages, arising out of or in connection with any errors in professional judgment or practice in the use of Nobel Biocare products. The clinician is also obligated to study the latest developments in regard to the Nobel Biocare products and their applications regularly. In cases of doubt, the clinician has to contact Nobel Biocare. Since the processing and application of this product are under the control of the clinician, they are his/her responsibility. We assume no liability whatsoever for damage arising therefrom.

Failure to observe the indicated limitations of use and working steps may result in failure.

Description: An individualized crown, cemented on natural teeth or abutments.

Materials: Lithium disilicate glass-ceramic

Indications:
 – Full contoured cement retained crown.
 – Applicable for anterior or posterior use, for functional and esthetic replacement of teeth.
 – Cemented on natural teeth or abutments.
 – Minimum crown thickness is indication dependant and the constraints incorporated within NobelProcera software (please refer to Ivoclar Vivident's IFU for additional information).

Contraindications:
 – Parafunctional tendencies, for example bruxisms and/or clenching.
 – Very deep sub gingival preparations.
 – Patients with substantially reduced residual dentition.

Procedural Precautions: Close cooperation between restorative dentist and dental technician is essential for a successful treatment. It is especially important to achieve proper stress distribution through adaptation and fitting of the crown or by adjusting the occlusion to the opposing jaw. In addition avoid excessive transverse loading forces, particularly in immediate loading cases. Because of the small size of prosthetic components, care must be taken that they are not swallowed or aspirated by the patient. Make sure that the minimum thickness indicated in the Ivoclar Vivident instructions for use are maintained even after polishing. If a try-in on the patient is carried out, we recommend cleaning the restoration under running water or with alcohol.

Sterility: The crown is delivered non-sterile and should be disinfected with a dedicated disinfection solution (e.g. 70% ethanol) before try-in and after try-in. After etching and conditioning of the crown, no additional solutions (e.g. disinfection solution) should be used.

Warnings:
 – Hexane is highly flammable and detrimental to health. Avoid contact of the material with skin and eyes. Do not inhale vapors and keep away from sources of ignition.
 – Do not inhale ceramic dust during finishing – use exhaust air discharge and mouth protection.
 – IPS Ceramic Etching Gel contains hydrofluoric acid. Contact with skin, eyes and clothing must be prevented at all costs, since the material is extremely toxic and corrosive. The etching gel is intended for professional use only and must not be applied introrsally (inside the mouth).

Handling procedures

Clinical Procedure

- Preparation of tooth**
 - Eliminate sharp edges, undercuts, and grooves.
 - Provide sufficient reduction (1.0 mm to 1.5 mm).
 - Provide adequate space for the crown.
 - Avoid sharp angles on the occlusal surface.
 - Avoid creating a deep fossa/cavity.
 - Provide sufficient (1.5 mm) occlusal/incisal reduction.
 - Avoid preparations that are excessively tapered or too close to parallel. The ideal total occlusal convergence is 6°-10°.
- Impression**
 - Take a conventional crown impression using your preferred impression material and method.
 - Send the impression to the laboratory.
- Shade determination**

Communicate both the shade of the preparation and the desired shade of the completed restoration to the dental technician. In this way, subsequent adjustments of the shade during cementation can be avoided.

Procedure to be carried out in the dental practice following delivery of the restoration from the dental laboratory.
- Conditioning of the preparation:**

Thoroughly clean the preparation once the temporary has been removed. Before it is conditioned, the restoration is tried in and the occlusion and articulation checked. If adjustments are required, the restoration may be polished exclusively in the restorative laboratory.

Condition the preparation in accordance with the selected cementation material.
- Conditioning of the restoration:**
 - Rinse restoration with water and blow dry with the air syringe.
 - Important:** Glass-ceramic materials must not be sandblasted!
 - Pre-treatment of the crown: Each inner surface with the hydrofluoric acid IPS Ceramic Etching Gel for 20 seconds, thoroughly with water and blow dry with the air syringe.
 - Apply the silane Monobond-plus to all inner surfaces for 60 seconds and air-dry. For conventional cementation, only etching is required, conditioning is not needed.
- Seating of the restoration**

Possibilities for esthetic cementation are decisive for the harmonious shade effect of an all-ceramic restoration. Depending of the indication, IPS e.max CAD restorations can be seated using either adhesive, self-adhesive or conventional cementation.

 - For the adhesive cementation of IPS e.max CAD restorations, Variolink® II, Variolink® Veneer or Multilink® Automix are the ideal cements.
 - For the self-adhesive cementation of IPS e.max CAD restorations, SpeedCEM is available.
 - It is not advisable to use classical/opaque phosphate (e.g. diamond-coated ceramic polishes – see light transmission of the all-ceramic and therefore compromise the esthetic appearance of allceramic restorations.*

Follow the instructions from the cementation supplier.
- Finishing and polishing:**
 - Occlusal adjustments are performed using fine diamonds (30 micron).
 - Polishing is carried out with ceramic polishing sets (e.g. diamond-coated ceramic polishes – see Ivoclar Vivident Chairside grinding recommendations for the dental office).
 - After removing all excess, fluorideize the tooth with e.g. Fluor Protector.
- Care notes**

As with natural teeth, precision milled IPS e.max CAD restorations require regular professional care. This is beneficial for both the health of the gingiva and teeth, as well as the overall appearance: The pumice-free Proxyl pink polishing paste is used to care for the surfaces without causing any wear. The low RDA value = 7 (Relative Dentin Abrasion) is a reliable confirmation to use a cleaning paste that is only a little abrasive. Scientific investigations and long-term clinical experience have proved the gentle effect compared to other pastes.

Laboratory Procedure

- Production:**
 - A model is made from the impression. Using the NobelProcera System, the model is scanned, the NobelProcera crown designed supported by supported by the tutorial found within the software. Once designed the order is sent to centralized production.
 - To achieve the optimal shade of the restoration, the blocks are selected on the basis of the desired tooth shade (Bleach BL or A-D), the determined die shade (ND 1-ND 9), and/or abutment shade, as well as the desired processing technique.
- Finishing:**
 - If adjustments are required, grind IPS e.max CAD restorations while in pre-crystallized (blue) state, if possible.
 - Only use suitable grinding instruments (see Ivoclar Vivadent grinding recommendations for labs), low rpms and light pressure to prevent delamination and chipping at the edges in particular.
 - Overheating of the glass-ceramic must be avoided.
 - The restoration is tried in on the dies and carefully finished.
 - Check proximal and occlusal contact points.
 - Make sure that the minimum thicknesses are maintained even after the minor adjustments.
 - Design surface textures.
 - Always clean the restoration with ultrasound in a water bath or blast with the steam jet before crystallization.
 - Make sure to thoroughly clean the restoration before further processing.
 - Do not blast the restorations with Al₂O₃ or glass polishing beads.
- Crystallization and Stain/Glaze Firing**

Refer to and follow the instructions provided by Ivoclar Vivadent for the IPS e.max CAD Crystallization and Stain/Glaze Firing procedures. Precision milled IPS e.max CAD crowns by NobelProcera can be finished in three different ways.

Option A: Crystallization and Glaze firing with IPS e.max CAD Crystall./Glaze paste
 In this processing technique, crystallization and Glaze firing are performed in one step. Characterizations is applied using IPS e.max CAD Crystall./Shades and Stains.

Option B: Crystallization and Stain/Glaze firing with IPS e.max CAD Crystall./Glaze spray
 IPS e.max CAD Crystall./Glaze Spray is used instead of IPS e.max CAD Crystall./Glaze Paste. The spray should only be used if IPS Object Fix auxiliary firing paste is adapted up to the margins.

Option C: Crystallization and separate Stain/Glaze firing using either IPS e.max CAD Crystall./Shades, Stains, Glaze or using IPS e.max Ceram Shades, Essences and Glaze
 In this processing technique, the IPS e.max CAD restorations are crystallized in a first step without applying any Stains and Glaze materials. Subsequently, the Stain and Glaze firing of the tooth-colored restoration is conducted.

Manufacturer: Nobel Biocare AB, Box 5190, SE-402 26 Göteborg, Sweden. Phone: +46 31 81 88 00. Fax: +46 31 16 31 52. www.nobelbiocare.com

Caution: The caution text "Federal (USA) law restricts the sale of this device to, or on the order of, a licensed physician or dentist" is shown on labels with "Rx Only".

Do not re-use. Re-use could result in local infection or transmission of infectious diseases

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