







INSTRUCTIONS FOR USE

71289 ML West 1102 Printed in Sweden. Date of issue 2011-02-04.
All rights reserved. Nobel Biocare, the Nobel Biocare logo and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Nobel Biocare. Product images in this folder are not necessarily to scale.

Connection	Platform	Titanium	Zirconia	Titanium	Zirconia
NobelReplace™ (internal tri-channel)	NP	X	X	X	X (Imp)
	RP	X	X	X	X (Imp)
	WP	X	X	X	X (Imp)
Brånemark System™, NobelSpeedy™ (external hex)	NP	X	X	X	X (Imp)
	RP	X	X	X	X (Imp)
	WP	X	X	X	X (Imp)
NobelActive™ (internal conical)	NP	X	X	X	X (Imp)
	RP	X	X	X	X (Imp)
	WP	X	X	X	X (Imp)
Straumann® Bone Level	NC	X (Ab)	X (Ab)	–	–
	RC	X (Ab)	X (Ab)	–	–
	4.8	X (Ab)	X (Ab)	X (Imp*)	–
Straumann® Octagon	6.5	X (Ab)	X (Ab)	X (Imp*)	–
	Yellow	X (Ab)	X (Ab)	X (Ab)	–
	Aqua	X (Ab)	X (Ab)	X (Ab)	–
AstraTech®	Lilac	X (Ab)	X (Ab)	X (Ab)	–
	3.3	X (Ab)	X (Ab)	X (Ab)	–
	3.8	X (Ab)	X (Ab)	X (Ab)	–
Camlog® Implant System	4.3	X (Ab)	X (Ab)	X (Ab)	–
	5.0/6.0	X (Ab)	X (Ab)	X (Ab)	–
	3.5	X (Ab)	X (Ab)	X (Ab)	–
Ankylos® Implant System	4.5	X (Ab)	X (Ab)	X (Ab)	–
	5.5	X (Ab)	X (Ab)	X (Ab)	–
	7.0	X (Ab)	X (Ab)	X (Ab)	–
3i®, Lifecore® Biomedical Restore, Sterngold® Implamed, Zimmer® Dental Taper-Lock	3.75	X	X	X	X (Imp)
	4.0	X	X	X	X (Imp)
	4.0	X	X	X	X (Imp)

(Ab) = Abutment level only (Imp) = Implant level only
* Implant level is supported through centralised scanning

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Disclaimer of liability: These prosthetic components are part of an overall concept and may only be used according to the instructions and recommendations of Nobel Biocare. Non-recommended use of products made by third parties in conjunction with Nobel Biocare prosthetic components will void any warranty or other obligation, express or implied, of Nobel Biocare. The user of Nobel Biocare prosthetic components has the duty to determine whether or not any product is suitable for the particular patient and circumstances. Nobel Biocare disclaims any liability, express or implied, and shall have no responsibility for any direct, indirect, punitive or other damages, arising out of or in connection with any errors in professional judgment or practice in the use of Nobel Biocare products. The clinician is also obliged to study the latest developments in regard to the Nobel Biocare products and their applications regularly. In cases of doubt, the clinician has to contact Nobel Biocare. Since the processing and surgical application of this product are under the control of the clinician, they are his/her responsibility. We assume no liability whatsoever for damage arising therefrom.

Description: A personalized bridge framework which after veneering at a dental laboratory, will be connected to endosseous dental implants or abutments intended for use as an aid in prosthetic rehabilitation.

Indications: NobelProcera™ and Procera™ Implant Bridge, is intended for multiple unit screw retained restorations at implant or abutment level. The NobelProcera™ and Procera™ Implant Bridge Zirconia should be designed to fit into a disc of Ø 60 mm x height 20 mm (illustration A). The connector dimension of a multi-unit framework depends on the distance between the supporting teeth (see table below and illustrations under B). The design is supported by real-time feedback through warming function in the NobelProcera™ Software.

Type - any position	Maximum length a=Arc length (mm) e=Extension length (mm)	Minimum connector and cross section height x width (mm) / area (mm²)
Free hanging arc	a ≤ 21.0	4.0 x 2.5 / 8.0
Free hanging arc	21.0 < a ≤ 35.0	4.0 x 3.0 / 9.4
Extension	e ≤ 10.0	4.0 x 3.0 / 9.4

Contraindications NobelProcera™ and Procera™ Implant Bridge, Zirconia:
– Cases where the mesial/distal cantilevers have a length of more than one unit.
– Cases with lengths that exceed the maximum limits.
– Bruxism.
– External radius less than 0.6 mm.

Instruction for clinician: We strongly recommend that clinicians, new as well as experienced implant users, always go through special training before undertaking a new treatment method. Nobel Biocare offers a wide range of courses for various levels of knowledge. For more info visit www.nobelbiocare.com. Always work with an experienced colleague the first time you employ a new treatment method. Nobel Biocare has a global network of mentors for this purpose.

Procedural Precautions: Close cooperation between surgeon, restorative dentist and dental laboratory technician is essential for a successful implant treatment. It is especially important to achieve proper stress distribution through passive adaptation and fitting of the bridge to the implant abutments; adjusting occlusion to the opposing jaw, avoiding excessive transverse loading forces, particularly in immediate loading cases. If the prosthesis metal substructure is made of gold alloy, this should have a high gold content. Because of the small size of prosthetic components, care must be taken that they are not swallowed or aspirated by the patient.

Sterility: The device is delivered non-sterile and should be cleaned and sterilized (if indicated) according to generally accepted procedures for dental laboratory work. Steam sterilize bridges and the prosthetic kit, if applicable for 5 min at 135° C/274° F.

Handling procedures

Clinical Procedure

- Place the impression coping onto the implant or Multi-unit Abutment and take an impression (C).
- Connect the healing abutment or temporary restoration.

Laboratory Procedure

- Produce a working model with removable gingival material (D).
- Create an acrylic frame work and a silicone index (E).
– A minimum material thickness of 0.9mm from the seating surface and 3mm up will be automatically added if needed (F).
- Using NobelProcera™ or Procera™ System – scan the model and acrylic framework. Design and order the Implant Bridge framework by the tutorial found within the software.

- After receiving the framework from Nobel Biocare production facility.**
- Finishing procedure for **Titanium** framework:
– Grind the titanium framework with hard metal burs only in one direction.
– Sandblast using 2.5bars of pressure utilizing 250µm aluminum oxide, at a 45° angle.
– Clean in an ultrasonic unit with distilled water for 5 minutes.
– Let the framework dry in room temperature 10-30 minutes before applying the bonding material.
– Apply titanium dental ceramic, acrylic or composite.
Note: Abutment screws/Prosthetic screws are not included. See screw chart for correct screw (table II and III).
 - Finishing procedure for **Zirconia** framework:
– If necessary, make minor adjustments using diamond impregnated finishing tools with fine grit size, under low pressure and using copious water irrigation.
– Sandblast using one bar of pressure utilizing 110µm aluminum oxide, at a distance of approximately 10mm.
– Apply zirconia dental ceramic.
Note: Abutment screws/Prosthetic screws are not included. See screw chart for correct screw (table II and III).

Clinical Procedure

- Connect the restoration onto the implant or abutment.
- Tighten the abutment (clinical) screws according to table I and the prosthetic Multi-unit screws 15Ncm, for Nobel Biocare systems use Screwdriver Unigrip™ and Manual Torque Wrench prosthetic and for other systems use appropriate screw driver and Torque Wrench (G and H).
Note: When using the NobelProcera™ or Procera™ Implant Bridge with external connection or compatible hex implants, use screw intended for the Brånemark System™ platform.

For additional information on restorative and dental laboratory procedures please consult treatment guidelines available at www.nobelbiocare.com or request latest printed version from a Nobel Biocare representative.

Materials: NobelProcera™ and Procera™ Implant Bridge Titanium – CP Titanium or Titanium alloy 90% Ti, 6% Al, 4% V.
NobelProcera™ Implant Bridge Zirconia – Yttria-stabilized Zirconiumoxide.
Abutment screws – Titanium alloy 90% Ti, 6% Al, 4% V.

Manufacturer: Nobel Biocare AB, Box 5190, SE-402 26 Göteborg, Sweden. Phone: +46 31 81 88 00. Fax: +46 31 16 31 52. www.nobelbiocare.com

Caution: The caution text "Federal (USA) law restricts the sale of this device to, or on the order of, a licensed physician or dentist" is shown on labels with "RX Only".

Straumann® is a registered trademark of Institut Straumann AG.
Camlog® is a trademark of Camlog Biotechnologies Group.
Astra Tech® is a trademark of Astra Tech Group.
3i® is a trademark of Biomet Group.
Ankylos® is a trademark of Dentsply Friadent Group.

Do not re-use. Re-use could result in local infection or transmission of infectious diseases.

III Abutment level

Multi-Unit Abutment prosthetic screw chart

	Prosthetic screw	Lab screw
NobelReplace™ NP	29285	29287
NobelReplace™ RP	29285	29287
NobelReplace™ WP, 6.0	29285	29287
Brånemark NP	29285	29287
Brånemark RP	29285	29287
Brånemark WP	29285	31163
NobelActive™ NP	29285	29287
NobelActive™ RP	29285	29287
Astra Tech™ Yellow	29285	29287
Astra Tech™ Aqua	29285	29287
Astra Tech™ Lilac	29285	29287
Astra Tech™ 4.5ST and 5.0ST	29285	29287
Straumann® Bone Level RC 3.3	29285	29287
Straumann® Bone Level RC 4.1/4.8	29285	29287
Straumann® Octagon soft tissue level 4.8	29285	29287
Straumann® Octagon soft tissue level 6.5	29286	31163
CAMLOG® 3.3 and 3.8	29285	29287
CAMLOG® 4.3	29285	29287
CAMLOG® 5.0 and 6.0	29285	29287

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