
NobelProcera™ Crown and Bridge

Crown and Bridge in Base Metal Alloy CoCr.

INSTRUCTIONS FOR USE

71281 ML West 1003 Printed in Sweden. Date of issue 2010-03-15.
All rights reserved. Nobel Biocare, the Nobel Biocare logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Nobel Biocare.
VITA VM13® is a registered trademark of VITA Zahnfabrik H. Rauter GmbH & Co. KG. Product images in this folder are not necessarily to scale.

Contraindications:
– Parafunctional tendencies, for example bruxisms and/or clenching
– If the coping or bridge is not veneered, and there is occlusal contact with other metal constructions/ restorations
– Cases with lengths that exceed the maximum limits
– Bridges must be designed to fit into a cylindrical disc with a diameter of 98mm x 15 (length x height)
– More than one free hanging pontic as cantilever

Procedural Precautions: Close cooperation between restorative dentist and dental laboratory technician is essential for a successful treatment. It is especially important to achieve proper stress distribution through adaptation and fitting of the crown by adjusting occlusion to the opposing jaw. In addition avoid excessive transverse loading forces, particularly in immediate loading cases. Because of the small size of prosthetic components, care must be taken that they are not swallowed or aspirated by the patient. If the restoration requires modification, pay attention to the marginal area and occlusal surfaces; a minimum thickness of 0.4mm should be kept. If modifying the restoration, use appropriate protection equipment to avoid inhale the Base Metal Alloy CoCr dust.
If a try-in on the patient is carried out, we recommend cleaning the restoration under running water or with alcohol.

Sterility: The coping or bridge is delivered non-sterile and should be disinfected with a dedicated disinfection solution (e.g. 70% ethanol) before try-in and again after try-in directly before cementation.

Handling procedures

Clinical Procedure

- Preparation of tooth:**
 - Eliminate sharp edges, undercuts, and grooves.
 - Provide sufficient reduction (1,2mm to 1,5mm).
 - Provide adequate space for the crown.
 - Avoid sharp angles on the occlusal surface.
 - Avoid creating a deep fossa/cavity.
 - Provide sufficient (1,5mm to 2mm) occlusal/incisal reduction.
 - Avoid preparation that are excessively tapered or too close to parallel. The ideal total occlusal convergence is 6–10°.
- Impression:**
 - Make an ordinary crown and bridge impression using your preferred impression material and method.
 - Send the impression to the laboratory.
Procedure after the restoration has been delivered from the dental laboratory
- Cementation:** Gently seat the restoration on the tooth and check both the occlusion and the interproximal contacts. The restoration should be in light occlusion. Excursive contact should be minimal.
- Removal:** To avoid overheating cross-cutting burs with copious irrigation during preparation should be used, if rotating instruments are used for removal.

Laboratory Procedure

- Production:** A model is made from the impression and scanned for producing a NobelProcera Crown or bridge with the NobelProcera Scanner and the NobelProcera software according to the tutorial found within the software. The framework should be designed with optimum veneering material support. The coping or bridge is made in a NobelProcera production site.
- Veneering:** After production, the crown is veneered with veneering material designed for use with CoCr with a CTE of 14.1. A veneering option is VITA VM13®.

Manufacturer: Nobel Biocare AB, Box 5190, SE-402 26 Göteborg, Sweden.
Phone: +46 31 81 88 00. Fax: +46 31 16 31 52. www.nobelbiocare.com

Caution: The caution text "Federal (USA) law restricts the sale of this device to, or on the order of, a licensed physician or dentist" is shown on labels with "Rx Only".

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ENGLISH

Important: Please read.


Disclaimer of liability: This prosthetic component is part of an overall concept and should only be used in conjunction with the recommended products according to the instructions and recommendation of Nobel Biocare. Use of non-recommended products made by third parties in conjunction with Nobel Biocare prosthetic components will void any warranty or other obligation, express or implied, of Nobel Biocare. The user of Nobel Biocare prosthetic components has the duty to determine whether or not any product is suitable for the particular patient and circumstances. Nobel Biocare disclaims any liability, express or implied, and shall have no responsibility for any direct, indirect, punitive or other damages, arising out of or in connection with any errors in professional judgment or practice in the use of Nobel Biocare products. The clinician is also obliged to study the latest developments in regard to the Nobel Biocare products and their applications regularly. In cases of doubt, the clinician has to contact Nobel Biocare. Since the processing and application of this product are under the control of the clinician, they are his/her responsibility. We assume no liability whatsoever for damage arising thereof. Failure to observe the indicated limitations of use and working steps may result in failure.

Description: An individualized coping or bridge cement retained solution on natural teeth or abutments, for use as the core structure for a crown.

Materials: Base Metal Alloy CoCr. Nickel, Cadmium and Beryllium free material. Material according to ISO 22674:2006.

Indications:
– Anterior or posterior coping or bridge for functional and esthetic replacement of teeth
– Tooth- and implant-supported restorations
– For a bridge span of up to 6 units
– Minimum thickness 0.4mm
– Any position in the mouth
– The connector dimension of a multi-unit framework depends on the distance between the supporting teeth (see table below). The design is supported by real-time feedback through warning function in the NobelProcera Software

Type – any position	Maximum length a=Arc length [mm] e=Extension length [mm]	Minimum connector and cross section height x width [mm] / area [mm²]
Free hanging arc	a ≤ 21,0 e ≤ 10,0	3,0 x 2,5 / 6,0 4,0 x 3,0 / 9,4

Contraindications:
– Parafunctional tendencies, for example bruxisms and/or clenching
– If the coping or bridge is not veneered, and there is occlusal contact with other metal constructions/ restorations
– Cases with lengths that exceed the maximum limits
– Bridges must be designed to fit into a cylindrical disc with a diameter of 98mm x 15 (length x height)
– More than one free hanging pontic as cantilever

Procedural Precautions: Close cooperation between restorative dentist and dental laboratory technician is essential for a successful treatment. It is especially important to achieve proper stress distribution through adaptation and fitting of the crown by adjusting occlusion to the opposing jaw. In addition avoid excessive transverse loading forces, particularly in immediate loading cases. Because of the small size of prosthetic components, care must be taken that they are not swallowed or aspirated by the patient. If the restoration requires modification, pay attention to the marginal area and occlusal surfaces; a minimum thickness of 0.4mm should be kept. If modifying the restoration, use appropriate protection equipment to avoid inhale the Base Metal Alloy CoCr dust.
If a try-in on the patient is carried out, we recommend cleaning the restoration under running water or with alcohol.

Sterility: The coping or bridge is delivered non-sterile and should be disinfected with a dedicated disinfection solution (e.g. 70% ethanol) before try-in and again after try-in directly before cementation.

Handling procedures

Clinical Procedure

- Preparation of tooth:**
 - Eliminate sharp edges, undercuts, and grooves.
 - Provide sufficient reduction (1,2mm to 1,5mm).
 - Provide adequate space for the crown.
 - Avoid sharp angles on the occlusal surface.
 - Avoid creating a deep fossa/cavity.
 - Provide sufficient (1,5mm to 2mm) occlusal/incisal reduction.
 - Avoid preparation that are excessively tapered or too close to parallel. The ideal total occlusal convergence is 6–10°.
- Impression:**
 - Make an ordinary crown and bridge impression using your preferred impression material and method.
 - Send the impression to the laboratory.
Procedure after the restoration has been delivered from the dental laboratory
- Cementation:** Gently seat the restoration on the tooth and check both the occlusion and the interproximal contacts. The restoration should be in light occlusion. Excursive contact should be minimal.
- Removal:** To avoid overheating cross-cutting burs with copious irrigation during preparation should be used, if rotating instruments are used for removal.

Laboratory Procedure

- Production:** A model is made from the impression and scanned for producing a NobelProcera Crown or bridge with the NobelProcera Scanner and the NobelProcera software according to the tutorial found within the software. The framework should be designed with optimum veneering material support. The coping or bridge is made in a NobelProcera production site.
- Veneering:** After production, the crown is veneered with veneering material designed for use with CoCr with a CTE of 14.1. A veneering option is VITA VM13®.

Manufacturer: Nobel Biocare AB, Box 5190, SE-402 26 Göteborg, Sweden.
Phone: +46 31 81 88 00. Fax: +46 31 16 31 52. www.nobelbiocare.com

Caution: The caution text "Federal (USA) law restricts the sale of this device to, or on the order of, a licensed physician or dentist" is shown on labels with "Rx Only".

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