**Gebrauchsanleitung: Straumann® CARES® Verschraubte Brücke und Straumann® CARES® Dolder® Steg**  
**Instructions for use: Straumann® CARES® Screw-retained bridge and Straumann® CARES® Dolder® Bar**  
**Mode d'emploi : Straumann® CARES® Bridge vissé et Straumann® CARES® Barre Dolder®**  
**Istruzioni per l'uso: Straumann® CARES® Ponte avvitato e Straumann® CARES® Barra Dolder®**

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**2. Placement in patient**  
 For additional information about the use of Straumann® products, call the Straumann customer service department. The following brochures provide additional guidance on the Straumann® Dental Implant System in relation to the bridge and bar devices.  
 – Crown and Bridge Restorations, Art. No. 152.255  
 – Retentive Systems for Implant-Borne Hybrid Dentures, Art. No. 152.252  
 – Information on the Surgical Procedures Art. No. 152.754  
 – Care and Maintenance of Surgical and Prosthetic Instruments Art. No. 152.008

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**English**

**Caution: U.S. Federal law restricts this device to sale by or the order of a licensed dentist.**

**Figure 1 – Straumann® CARES® Screw-retained bridge**

**Figure 2 – Straumann® CARES® Dolder® Bar**

**Product description**  
 The Straumann® CARES® Screw-retained bridge (bridge device) and Straumann® CARES® Dolder® Bar (bar device) devices are used for the restoration of Straumann® dental implants with different endosteal diameters, lengths and platforms.

The bridge and bar devices allow for individual customization regarding function and esthetics. They attach directly to a Straumann® dental implant, meaning that the abutment/secondary part is contained in the device. The devices are intended to be finished into a bridge or denture using standard dental laboratory techniques and materials. The bridge and bar devices are CAD-derived, CAM-produced and have a scanner as its data source. The devices are made entirely of cobalt-chromium alloy (coron®).

These instructions for use are valid for Straumann® CARES® Screw-retained bridge and Straumann® CARES® Dolder® Bar devices of the following prosthetic lines: RN (Regular Neck) Ø 4.8 mm and

WN (Wide Neck) Ø 6.5 mm. The device is currently not available for Straumann® Soft Tissue level implants with a NN (Narrow Neck) Ø 3.3 mm prosthetic interface and Straumann Bone level implants with a NC (Narrow CrossFit®) Ø 3.3 mm and RC (Regular CrossFit®) Ø 4.1 mm prosthetic interface.

The device is delivered from Straumann to the receiving party with the following content:  
 – Bridge or bar device  
 – Final patient screws  
 – Patient labels  
 – Instructions for use

**Indications for use**  
 The bridge and bar devices are indicated to be used only with implants from the Straumann® Dental Implant System. The bridge and bar devices are indicated for screw-retained restorations on Straumann® implants with prosthetic interface RN (Regular Neck) Ø 4.8 mm and WN (Wide Neck) Ø 6.5 mm.  
 The bridge device is intended to be used as a bridge framework in the treatment of partially or totally edentulous jaws for the purpose of restoring chewing functions.  
 The bar device is intended to be used as a retentive element in the treatment of edentulous jaws in conjunction with an overdenture for the purpose of restoring chewing functions.

Consult the brochure "Basic Information on the Surgical Procedures" (Art. No. 152.754) for more information on indications of implants.

**Contraindications**  
 – Allergies or hypersensitivity to chemical ingredients of materials used: coron® (cobalt-chromium alloy, Ni and Be free alloy, type 4, according to ISO 22674) for bridge and bar devices, TAN (titanium alloy Ti-6Al-7Nb) for the screw  
 – Usage on any other implant line than the Straumann® Dental Implant System  
 – In case of implant divergences greater than 40° between two adjacent implants  
 – Consult the brochure "Basic Information on the Surgical Procedures" (Art. No. 152.754) for information on indications and contraindications of implants, such as required minimum number of implants, implant type, diameter, loading protocols.  
 – Patients with bruxism are not to be treated with a bridge and bar device as an overload of the device may occur.

**Warning**  
 The following descriptions are not sufficient for immediate use of the Straumann® Dental Implant System. Knowledge of dental implantology and instruction in the handling of the Straumann® Dental Implant System provided by an operator with the relevant experience are always necessary. Courses and training sessions for Straumann® Dental Implant System are regularly offered by the ITI (International Team for Implantology) and Institut Straumann AG.

**Precautions**  
 – All products must be secured against aspiration when used intraorally.  
 – In case of wrong usage implant soft, soft tissue irritation or injury to the patient may be the consequence.  
 – Always use the new patient screws for final device insertion in the patient that is included in your shipment.  
 – Ensure regular relining of the denture to avoid an overload of the bar device.  
 – Ensure regular cleaning of the denture to avoid inflammation of the soft tissue.

For further information on the use of the Straumann® CARES® Dolder® bar, it is recommended to combine the device only with original matrices of the Straumann® CARES® Dolder® system.  
 – When reducing the height of the bar coping, make sure to reduce it from the top to the level of the bar. A minimum height of the bar coping of 3 mm must be maintained (Figure).

**5. Placement in patient**  
 The bridge and bar devices are inserted into the implant without applying cement. They are inserted and tightened with the appropriate screw using a Straumann® SCS screwdriver, ratchet and torque control device. Check that the SCS screwdriver is not defect to ensure appropriate usage and aspiration safety of the screw.

| Abutment platform | Additional screw | Tightening force  |
|-------------------|------------------|---|
| RN<br>WN          | 048.356          | – on implant analogs:<br>handtight<br>– try-in in patient: handtight<br>– final placement: 35 Ncm |

**Cleaning and disinfection**  
 Cleaning should be done with distilled water in ultrasound or with ethyl alcohol. Do not immerse coron® frameworks into a pickling bath. Consult the brochure "Care and Maintenance of Surgical and Prosthetic Instruments" (Art. No. 152.008) for additional information.

**Sterilization**  
 Straumann® abutments and components are not sterile when delivered. Straumann recommends the following procedure for sterilization prior to use:

| Material | Method                | Conditions              |
|----------|-----------------------|-------------------------|
| coron®   | Autoclave, moist heat | 134 °C [273 °F], 18 min |
| TAN      | Autoclave, moist heat | 134 °C [273 °F], 18 min |

**Note**  
 Parts that have been modified or altered in any way may require a different sterilisation procedure.


**Instructions for use**  
 For detailed instructions consult "Basic information on Straumann® CARES® Dolder® Bar procedures" and "Basic information on Straumann® CARES® Screw-retained bridge procedures".

**1. Digitalization of patient situation**  
 Prequisite for a good fitting final device is an accurate impression and subsequent digitalization. Digitize the patient model using the Straumann® CARES® Scan CS2 scanner.

**2. CAD of Straumann® CARES® Screw-retained bridge and Straumann® CARES® Dolder® Bar**  
 Design the bridge and device in the Straumann® CARES® Visual software. Beware that inappropriate designs may cause failure of the device and/or implant harm to the patient. The software does not allow to fall below a minimum value for cross-section areas or to exceed the number of pontics (maximum 4) and freends (maximum 1). When designing bridges with freends (1) or bridges with the maximum number of pontics (4) in the front region, it is required to design the cross-section areas accordingly to ensure a bridge design with sufficient strength as the minimal cross-section defined by the software might not be sufficient in these extreme cases.

The Straumann Product is part of an overall concept and must be used only in conjunction with the corresponding original components and instruments distributed by Institut Straumann AG. Its ultimate parent company and all affiliates or subsidiaries of such parent company ("Straumann"), except if stated otherwise in these instructions for use. If use of products made by third parties is not recommended by Straumann in these instructions for use, any such use will void any warranty or other obligation, express or implied, of Straumann.

**Material-specific design recommendations**  
 To ensure a bridge framework with sufficient strength, take A, B, W, H and C into consideration during the design process.

**Figure:** Pontic situation [A], extensions situation [B], height [H], width [W] and cross-section area [C] of connector.

**Pontic situation A**  
 Figure 3 shows recommendations for height [H], width [W] and cross-section area [C] of connectors for a bridge that is placed from a tooth position X to a tooth position Y.

Move along the lines of tooth positions to determine dimensions H, W and C of the connector.  
 Example: For a pontic situation from tooth position 11 to position 13 it is recommended to design all connectors inbetween with a height of at least 3.0 mm, a width of at least 2.0 mm and cross-section area of at least 5.1 mm<sup>2</sup>.

**Extensions situation B**  
 Figure 4 shows recommendations for height [H], width [W] and cross-section area [C] of the connector for a bridge framework with an extension. The extension is between tooth position X and tooth position Y.

Move along the lines of tooth positions of the extension to determine dimensions H, W and C of the connector.  
 Example: For an extension situation with a connector between tooth position 11 and tooth position 12, it is recommended to design the connector with a height of at least 4.5 mm, a width of at least 3.0 mm and a cross-section area of at least 12.6 mm<sup>2</sup>.

**3. CAM of Straumann® CARES® Screw-retained bridge and Straumann® CARES® Dolder® Bar**  
 The device is produced in the Straumann milling center according to your CAD file.

**4. Finishing of device**  
**a. Unpacking** – Each package contains all components for one patient case. Take care when unpacking the device as it may have sharp edges and small parts may fall out.  
**b. Cleaning, try-in on model and patient** – Clean the device before placing it on implant analogs or implants. Try on the device on the master cast model and then in the patients mouth for accurate fit. Always use corresponding lab screws (order separately) until the final placement of the device. When shipping the device always ensure that it is screwed handtight on the model to avoid any harm of the prosthetic connection during transport.  
**c. Sandblasting/polishing** – Use polishing aids to protect the prosthetic connection as well as the screw and the screw channel when sandblasting and/or polishing the device. Use appropriate pressure for sandblasting to avoid reduction of the device cross-section. Use standard polishing methods/materials (110 µm Al<sub>2</sub>O<sub>3</sub>, 2 bar pressure).  
**d. Veneering** – Clean the device before veneering. All commercially available bonding ceramics that are intended for veneering of dental cobalt-chromium alloys and that are suitable for a CTE value of 14.4 x 10<sup>6</sup>/K can be used (such as HeroCeram®). The processing instructions of the veneering material manufacturers apply. Protect the prosthetic connection with a polishing aid as well as the screw and screw channel. The surface must be cleaned thoroughly by brushing with distilled water after each ceramic firing to remove soluble oxides.

**5. Placement in patient**  
 The bridge and bar devices are inserted into the implant without applying cement. They are inserted and tightened with the appropriate screw using a Straumann® SCS screwdriver, ratchet and torque control device. Check that the SCS screwdriver is not defect to ensure appropriate usage and aspiration safety of the screw.

| Abutment platform | Additional screw | Tightening force  |
|-------------------|------------------|---|
| RN<br>WN          | 048.356          | – on implant analogs:<br>handtight<br>– try-in in patient: handtight<br>– final placement: 35 Ncm |

**Cleaning and disinfection**  
 Cleaning should be done with distilled water in ultrasound or with ethyl alcohol. Do not immerse coron® frameworks into a pickling bath. Consult the brochure "Care and Maintenance of Surgical and Prosthetic Instruments" (Art. No. 152.008) for additional information.

**Sterilization**  
 Straumann® abutments and components are not sterile when delivered. Straumann recommends the following procedure for sterilization prior to use:

| Material | Method                | Conditions              |
|----------|-----------------------|-------------------------|
| coron®   | Autoclave, moist heat | 134 °C [273 °F], 18 min |
| TAN      | Autoclave, moist heat | 134 °C [273 °F], 18 min |

**Note**  
 Parts that have been modified or altered in any way may require a different sterilisation procedure.

**Instructions for use**  
 For detailed instructions consult "Basic information on Straumann® CARES® Dolder® Bar procedures" and "Basic information on Straumann® CARES® Screw-retained bridge procedures".

**1. Digitalization of patient situation**  
 Prequisite for a good fitting final device is an accurate impression and subsequent digitalization. Digitize the patient model using the Straumann® CARES® Scan CS2 scanner.

**2. CAD of Straumann® CARES® Screw-retained bridge and Straumann® CARES® Dolder® Bar**  
 Design the bridge and device in the Straumann® CARES® Visual software. Beware that inappropriate designs may cause failure of the device and/or implant harm to the patient. The software does not allow to fall below a minimum value for cross-section areas or to exceed the number of pontics (maximum 4) and freends (maximum 1). When designing bridges with freends (1) or bridges with the maximum number of pontics (4) in the front region, it is required to design the cross-section areas accordingly to ensure a bridge design with sufficient strength as the minimal cross-section defined by the software might not be sufficient in these extreme cases.

The Straumann Product is part of an overall concept and must be used only in conjunction with the corresponding original components and instruments distributed by Institut Straumann AG. Its ultimate parent company and all affiliates or subsidiaries of such parent company ("Straumann"), except if stated otherwise in these instructions for use. If use of products made by third parties is not recommended by Straumann in these instructions for use, any such use will void any warranty or other obligation, express or implied, of Straumann.

**Material-specific design recommendations**  
 To ensure a bridge framework with sufficient strength, take A, B, W, H and C into consideration during the design process.

**Validity**  
 Upon publication of these instructions for use, all previous versions are superseded.

© Institut Straumann AG, 2010. All rights reserved. Straumann® and/or other trademarks and logos from Straumann® mentioned herein are the trademarks or registered trademarks of Straumann Holding AG and/or its affiliates.

Dolder® is a trademark registered by Prof. Dolder, exclusively licensed to Cendres+Métaux Holding SA.

**Français**

**Figure 1 – Straumann® CARES® Bridge vissé**

**Figure 2 – Straumann® CARES® Barre Dolder®**

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